

**ARGON MEDICAL DEVICES, INC.**  
1445 Flat Creek Road  
Athens, Texas 75751 USA  
Tel: +1 800-927-4669;  
Tel: +1 903-675-9321  
www.argonmedical.com

**EMERGO EUROPE**  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
+31 70 345 8570

**Argon Medical Devices UK Ltd**  
Eastgate Business Centre  
Eastern Avenue  
Burton-on-Trent  
Staffordshire  
DE13 0AT

<https://www.argonmedical.com/resources/product-information>

X9585354/0321A

English	The symbols glossary is located electronically at <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Bulgarian	Речникът на символите е наличен в електронен вид на адрес <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Czech	Slovník symbolů v elektronické formě je umístěn na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Danish	Symbolforklaringen findes elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
German	Das Glossar der Symbole kann unter <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Spanish	El glosario de símbolos puede consultarse en línea en la dirección <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Finnish	Elektroninen symbolisanasto on osoitteessa <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
French	Le glossaire des symboles est disponible au format électronique à l'adresse <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Latvian	Simbolių žodyno elektroninė forma pateikta <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Lithuanian	Simbolų glosarijs ir atrodams elektroniskā veidā vietnē <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Norwegian	Symbolforklaringen befinner seg elektronisk på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Romanian	Glosarul simbolurilor este stocat electronic la <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Russian	Электронный глоссарий символов расположен по адресу <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Slovenian	Slovar simbolov se nahaja v elektronski obliki na <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> .
Swedish	Symbolforklaringen hittas elektroniskt på <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a>
Turkish	Semboller sözlüğü elektronik ortamda <a href="http://www.argonmedical.com/symbols">www.argonmedical.com/symbols</a> adresinde yer almaktadır

## ENGLISH

### Bone Marrow Harvest Needle

**Intended Use/Purpose:** The Bone Marrow Harvest Needle are designed for allogenic and autologous bone marrow harvest-transplantation procedures.

**Devices Description:** The Bone Marrow Harvest Needle offers a standard bevel tip or a double-diamond tip and side holes to allow for efficient aspiration. The product features an ergonomic designed twist-lock handle provides greater clinical comfort and luer lock connector that allows for secure syringe attachment.

**Indications for Use:** Intended for the purpose of Harvesting Bone Marrow.

**Contraindications:** For use only for harvesting bone marrow. These needles should be used only by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of harvesting bone marrow. Physician judgement is required when considering harvest on patients with bleeding disorder or receiving anti-coagulant medications.

**Potential Complications/Risks :** It is not uncommon for patients to experience the following :

- Nerve injury including puncture of the spinal cord or nerve roots potentially resulting in radiculopathy, paresis or paralysis.
- Thrombus or other materials resulting in symptomatic pulmonary embolism or other clinical sequelae.
- Hemothorax or pneumothorax.
- Infection.
- Unintended puncture wounds including vascular puncture and dural tear.
- Bleeding or hemorrhage
- Hematoma
- Pain
- Discomfort
- Allergic reaction

Pulmonary and cardiovascular related complications occur mostly in rare situations when the device is used for sternal aspirations.

When use as intended by a physician, the benefits of using the Bone Marrow Biopsy Needle outweigh the risks associated with the use of the device.

#### Warning

- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the component are damaged during procedure.
- Do not bend the needle excessively, it may result damage/breakage to the device.
- This device was designed, tested and manufactured for single patient use only.
- When used repeatedly in the same patient, the device should be inspected for damage or wear after each core sample is taken. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness and/or injury. DO NOT reuse, reprocess or re-sterilize this device.
- Reconditioning, refurbishing, repair, modification, of the device to enable further use is expressly prohibited.

#### Precaution

- This device should be manipulated only while under fluoroscopic observation with radiographic equipment that provides high quality images.
- Access to the vertebral body via pedicle requires a minimum pedicle width of 5mm.
- Knowledge of the dimensions of the site of insertion as assessed by MRI, CT, or other imaging method.

**Note:** These instructions are not meant to define or suggest any medical or surgical technique. The individual physician is responsible for the proper procedure and techniques to be used with this device.

#### Harvest Procedure:

1. After suitable anesthesia is achieved, place the patient in the lateral right or left lateral position.
2. Using sterile technique, prepare the skin with antiseptic and drape.
3. Hold the needle with the proximal end in palm and the index finger against the shaft near the tip. This position stabilizes the needle and allows for better control.
4. Introduce the needle through the skin and bring it into contact with the posterior iliac crest.
5. Using gentle, but firm pressure, advance the needle, rotating it in an alternating clockwise/counterclockwise motion. Entrance into the marrow cavity is generally detected by decreased resistance. (All of the side holes at the distal end of the needle must be introduced into the marrow cavity beyond the cortical bone, otherwise air and extra bony soft tissue may appear with the aspirated marrow.)
6. Once the needle is in place, remove the stylet by rotating the upper section 90°, and pulling straight out.
7. Attach a syringe with a luer taper to the hub of the bone marrow harvest needle using a firm push and twist motion.
8. Apply suction by withdrawing the syringe plunger. Remove the syringe with the harvested marrow.

- Repeat the harvest procedure until an appropriate amount of marrow is obtained to satisfy the clinical requirement.

**Disposal:** After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

**Storage:** Store at controlled room temperature.

**NOTE:** In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at [quality.regulatory@argonmedical.com](mailto:quality.regulatory@argonmedical.com) as well as to the competent health authority where the user/patient resides.

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