ARGON MEDICAL DEVICES, INC.
 1445 Flat Creek Road
 Athens, Texas 75751 USA
 Tel: +1 (903) 675 9321
 Tel: +1 (800) 927 4669
 www.argonmedical.com

EMERGO EUROPE
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 +31 70 345 8570

Argon Medical Devices UK Ltd
 Eastgate Business Centre
 Eastern Avenue
 Burton-on-Trent
 Staffordshire
 DE13 0AT

RxOnly MD REF LOT

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovník symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Glossar der Symbole kann unter www.argonmedical.com/symbols abgerufen werden
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse en línea en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Latvian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Lithuanian	Simbolu glosarijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Russian	Электронный глоссарий символов расположен по адресу www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Arabic	www.argonmedical.com/symbols بتوفر مسرد الرموز إلكترونياً على

IFU6600-0221/0221A

Intended Use /Purpose

The Hystero-Salpingography Catheter (HSG) is designed for use in hysterosonography or hysterosalpingography procedures.

Device Description

The HSG Catheter balloon provides secure placement in the uterus or cervical canal and effectively seals to prevent retrograde leakage of saline or contrast media. The malleable introduction sheath retains curvature to facilitate easy placement into the cervical canal and it also include a fixed volume syringe to prevent from over-inflation of the balloon.

It has color-coded stopcocks on both catheter lumens with the following functions:

Blue Stopcock – is for balloon inflation.

White Stopcock – is for saline/contrast media infusion.

The HSG catheters are available in the following sizes with corresponding balloon diameters with air and saline.

Catalog No.	Accepts GW	Balloon Dia. W/ air	Balloon Dia. W/ saline
660005040 (5.5Fr)	0.035" (0.89mm)	11mm with 2.0ml air	10mm with 0.8ml saline
660007040 (7Fr)	0.038" (0.96mm)	13mm with 3.0ml air	12mm with 1.2ml saline

Indication for Use

HSG Catheters are intended for use in the injection of contrast material in the examination of the uterus and fallopian tubes.

Contraindications for Use

Physician judgment should be used in determining in which patients Hystero-Salpingography or Hystero-Sonography procedures are contraindicated. These procedures are not recommended for patients with active pelvic infections.

Side Effects

- It is not uncommon for patients to experience mild pain during or after the procedure, bleeding, exposure to radiation, pelvic infection, and inflammation of the endometrium or the fallopian tubes.

Duration

Transient, less than 60 minutes.

Warnings

- Do not use this device with a power injector for the delivery of the contrast media.
- This device was designed, tested and manufactured for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Do not continue to use if any of the components are damaged during the procedure.
- Do not use Lipiodol or other oil-based contrast media.
- Do not apply excessive force or deflection, this could result in damage and/or breakage of the device.

Risks

- Associated risks include perforation, fallopian tube infection or an allergic reaction due to contrast injection.

Precautions

- These devices should be used only by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of Hystero-Salpingography and Hystero-Sonography.
- A contrast and saline will be delivered by injection through the catheter as part of the procedure and the delivery rate is controlled by physician.

Pre-Procedure Preparation

- Fill the catheter with contrast liquid and close the stopcock. Remove the contrast syringe to facilitate handling of the catheter.
- Check that the blue stopcock on the balloon port is open and mount the syringe with the required quantity of air (X- ray) or sterile saline (ultrasound).
- If sterile saline is used in the balloon the air is deflated from the system by filling and emptying the balloon with sterile saline before use
- Inflate and deflate the balloon (extend and retract the plunger on the syringe) to insure its functionality and integrity.

Procedure

- Insert the introduction tube (IT) into the vagina and guide the balloon catheter through the (IT) into uterus (max 6 cm).
- Inflate the balloon and close the stopcock. Pull back the catheter to secure the balloon thus preventing outflow of contrast liquid into the vagina.
- Fix the introducer tube onto the reinforced part of the catheter shaft at the bifurcation.

- Mount the contrast syringe, open the stopcock, inject contrast liquid and carry out the examination.
- When the examination has been completed, deflate the balloon by opening the balloon port stopcock and drawing back the syringe. Remove the catheter if necessary using contrast injection and fluoroscopy in order to present the cervix uteri.

Disposal for Use

After use, this product may be potential biohazard. Dispose in accordance with hospital policies and procedures for waste and biohazard materials.

Storage

Store at standard ambient temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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