ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd

Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

IFU9012/0321A

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Ein elektronisches Symbolglossar finden Sie unter www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

ENGLISH

Needles

Intended Use/Purpose:

Argon Needles are designed for percutaneous punctures for guide wire introduction.

Devices Description:

The products addressed by these instructions for use includes : AMC/3 Needle, AMC/4 Needle, Courmand-Style Needle, Percutaneous Entry Needle, and Winged Arterial Needle.

This product offers either one piece, two piece or three piece assemblies, and are available in Needle size : 18ga, 19ga, or 21ga with length of 2.5cm, 3.8cm, 5.1cm, 5.7cm, 7cm, 7.3cm or 5.4cm.

Indication for Use:

Argon Needles are indicated for used in percutaneous entry procedure.

Warnings:

- Contents are supplied sterile and are intended for single use only. Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess, or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged or if the expiry date has been exceeded.
- Do not continue to use if any of the components are damaged during the procedure.

Caution: Remove Winged Needle Adaptor before attempting to pass a guide wire through the hub of the needle or before attaching a syringe.

The guidewire should not be withdrawn through the needle. Damage or shearing of guidewire may occur.

Precautions :

- This device is intended for use by trained, qualified licensed physician.

Direction for Use:

- Prepare and drape the skin at the intended puncture site; achieve local anesthesia as necessary.
- Using imaging guidance insert the needle through the skin
- Using imaging guidance advance the needle to the intended vessel or target site.
- Advance the appropriate guidewire of choice through the needle.
- Confirm access to desired location with wire
- Remove needle and continue with procedure.

Possible complications :

- Minor Pain at insertion site
- Minor Bleeding
- Hematoma and insertion site
- Infection
- Vessel puncture

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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