ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road, Athens, Texas 75751 USA
Tel: +1 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

Argon Medical Devices UK Ltd
Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

X9585765/Rev.0321A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Ein elektronisches Symbolglossar finden Sie unter www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusán itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

Intended Use/Purpose : Co-Axial Introducer Needle is used with either a side notch or full core biopsy needle device in performing soft tissue percutaneous biopsies.

Device Description : The Co-Axial Introducer Needle features an echogenic tip for accurate placement under ultrasound guidance, and a centimeter markings that provide reference for depth placement. The hubs are color coded for easy gauge size determination and it has an adjustable needle stop that allows the user to restrict the forward movement, effectively localizing the needle tip to the biopsy site.

Indications for Use : Co-Axial Introducer Needle is intended for a range of biopsy applications of soft tissue.

Warning

- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury. Do not reuse, reprocess or re-sterilize this device.
- Inspect the package integrity before use.
- Do not use if package is open or damaged or if the expiry date has been exceeded.
- The instrument should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of core needle biopsy.
- Use only for core biopsies of soft tissue as determined by a licensed physician. Attempting to biopsy hard tissue or masses may reduce the performance, causing damage to the instrument or potentially harm the patient. This instrument is intended for soft tissue biopsy only.
- Incorrect setting of stroke length on the device may cause a sample to be incorrectly taken.
- Handle in a manner, which will prevent accidental puncture.

Note: The device compatibility been tested with Argon Biopsy products. Other manufacturer's biopsy instrument has not been tested for compatibility.

Precaution

- Verify integrity of needle before use. If needle is damaged, appropriately discard the entire device (see disposal section) and prepare a new instrument.
- The introduction of the needle into the body should be carried out under imaging control (ultrasound)
- Do not bend the needle.
- Do not contact the needle tips to a hard surface while expelling a sample.

Possible/Side Effect

- It is not uncommon for patients to experience:
 - Bleeding, infection at insertion site
 - Fever
 - Pain at the biopsy site
 - Swelling at the biopsy site,
 - Adjacent tissue injury
 - Hematoma.

- ***When used as intended by a physician, the benefits of using the Biopsy Instrument Needle outweigh the risks associated with the use of the device.

DIRECTIONS FOR USE:

1. Prepare the biopsy site appropriately and locate the target area.
2. Advance a Co-axial Introducer Needle, using standard technique, to a position proximal to the lesion. Position the adjustable needle stop proximal to the skin.
3. Remove the stylet from the Co-axial Introducer Needle, leaving the cannula in position.
4. Choose your complementary biopsy needle.
5. Insert the biopsy needle through the coaxial needle.
6. Fire and retrieve the biopsy sample with the biopsy needle.
7. Removed the biopsy needle.

Disposal: After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

Storage: Store at controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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