ARGON

MEDICAL DEVICES

ARGON MEDICAL DEVICES, INC.

1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 800-927-4669;
Tel: +1 903-675-9321
www.argonmedical.com

Argon Medical Devices UK Ltd

Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE13 0AT

EMERGO EUROPE

Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

20-2007-52 /0321A

<https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgarian	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů je elektronicky umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Ein elektronisches Symbolglossar finden Sie unter www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en www.argonmedical.com/symbols
Estonian	Sümbolite seletuste kogu on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Elektroninen symbolisanasto on osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible au format électronique à l'adresse www.argonmedical.com/symbols
Croatian	Pojmovnik simbola nalazi se u elektroničkom obliku na www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodyno elektroninė forma pateikta www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen befinner seg elektronisk på www.argonmedical.com/symbols
Polish	Glosariusz symboli w postaci elektronicznej dostępny jest pod adresem www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente no endereço www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Slovenian	Slovar simbolov se nahaja v elektronski obliki na www.argonmedical.com/symbols .
Swedish	Symbolforklaringen hittas elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır

Intended Use/Purpose

Introducer sheath and PTC Catheter Needles are used in introduction for drainage procedures.

Devices Description

The Introducer Sheath Needle offers radiopaque sheaths for guidewire insertion or fluid aspiration. The needle has a trocar-style tip for easy insertion and it's available in 15cm and 20cm lengths.

The PTC Catheter Needles are available in 5F diameter in 17 and 25cm lengths. It consists of a polyethylene introducer sheath over a stainless steel trocar tip or bevel tip cannula.

Indication for Use

Introducer Sheath Needle is for the introduction of a .038" guidewire used in abscess, biliary, nephrostomy, and other fluid collection drainage procedures.

PTC Catheter Needles are intended for introducing a 0.038" guidewire, as a secondary puncture site for biliary drainage procedures, for contrast injection and pathfinding.

Warnings :

- Contents are supplied sterile and are intended for single use only. Do not re-sterilize, reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Inspect the package integrity before use.
- Do not use if package is open or damaged or if the expiry date has been exceeded.
- Do not continue to use if any of the components are damaged during the procedure.
- These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

Possible Complication: It is not uncommon for patients to experience:

- Injury at site
- Bleeding, infection
- Pain & hematoma.

Precautions

- This device is intended for use by qualified licensed physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of fluid collection.

Pre-Procedural Preparations :

- Prepare the Introducer Sheath Needle assembly by inserting the stylet needle into the introducer sheath and locking the luer assemblies. Ensure all connections are tight without overtightening. Overtightening may result in damage to components.

Direction for Use :

- Prepare and drape the skin at the intended puncture site; achieve local anesthesia as necessary.
- Make an adequate skin incision at the determined entry site.
- Advance the introducer sheath with stylet needle assembly through the skin.
- Using imaging guidance, advance the needle assembly to the desired location.
- Remove the stylet needle leaving the introducer sheath in place.
- Advance the appropriate guidewire of choice through the introducer sheath.
- Confirm access to desired location with wire.
- Remove sheath over wire and continue with procedure.

Storage: Store at controlled room temperature.

Disposal: After use, this product may be a potential biohazard. Dispose in accordance with applicable laws and regulations.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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