Screw Dilator

ARGON MEDICAL DEVICES, INC.
1445 Flat Creek Road
Athens, Texas 75751 USA
Tel: +1 (903) 675 9321
Tel: +1 (800) 927 4669
www.argonmedical.com

EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands
+31 70 345 8570

CE 2797

Argon Medical Devices UK Ltd
Eastgate Business Centre
Eastern Avenue
Burton-on-Trent
Staffordshire
DE 13 0AT

PN 20-2007-51/0321A

STERILE EO  MD RxOnly

 REF  LOT

 <https://www.argonmedical.com/resources/product-information>

English	The symbols glossary is located electronically at www.argonmedical.com/symbols
Bulgaria	Речникът на символите е наличен в електронен вид на адрес www.argonmedical.com/symbols
Czech	Slovníček symbolů v elektronické formě je umístěn na www.argonmedical.com/symbols
Danish	Symbolforklaringen findes elektronisk på www.argonmedical.com/symbols
German	Das Symbolglossar ist elektronisch verfügbar unter www.argonmedical.com/symbols
Greek	Το γλωσσάρι των συμβόλων παρέχεται ηλεκτρονικά στη διεύθυνση www.argonmedical.com/symbols
Spanish	El glosario de símbolos puede consultarse electrónicamente en la dirección www.argonmedical.com/symbols
Estonian	Sümbolite sõnastik on elektrooniliselt saadaval aadressil www.argonmedical.com/symbols
Finnish	Symbolisanasto on sähköisessä muodossa osoitteessa www.argonmedical.com/symbols .
French	Le glossaire des symboles est disponible sous forme électronique à l'adresse www.argonmedical.com/symbols
Hungarian	A szimbólumok összefoglalása elektronikusan itt található: www.argonmedical.com/symbols
Italian	Il glossario dei simboli è disponibile in formato elettronico all'indirizzo www.argonmedical.com/symbols
Lithuanian	Simbolių žodynas elektronine forma pateiktas www.argonmedical.com/symbols
Latvian	Simbolu glosārijs ir atrodams elektroniskā veidā vietnē www.argonmedical.com/symbols
Dutch	Uitleg over de symbolen is elektronisch beschikbaar op www.argonmedical.com/symbols
Norwegian	Symbolforklaringen finnes seg elektronisk på www.argonmedical.com/symbols
Polish	Wersja elektroniczna słownika symboli znajduje się na stronie www.argonmedical.com/symbols
Portuguese	O glossário de símbolos está localizado eletronicamente em www.argonmedical.com/symbols
Romanian	Glosarul simbolurilor este stocat electronic la www.argonmedical.com/symbols
Slovak	Glosár symbolov je v elektronickej podobe k dispozícii na stránkach www.argonmedical.com/symbols
Swedish	Symbolordlistan finns elektroniskt på www.argonmedical.com/symbols
Turkish	Semboller sözlüğü elektronik ortamda www.argonmedical.com/symbols adresinde yer almaktadır
Arabic	يوفر مصدر الرموز إلكترونياً على www.argonmedical.com/symbols

Intended Use/Purpose

The Screw Dilator product family is intended for tissue dilation.

Device Description

Screw dilators are intended for tissue dilation in non-vascular procedures, allowing dilation up to 14 French by using only one dilator. They feature smooth rounded tips with gradual tapering for facilitating tissue penetration, radiopaque for fluoroscopic visualization, and a luer lock hub for contrast injection and guidewire exchange.

Indication for Use

The Screw Dilator is used in introduction and placement of guidewires and catheters for biliary, nephrostomy, and general percutaneous fluid collection and drainage.

Contraindications for Use

There are no contraindications associated with the Screw Dilators. However, assessment of patient overall physical state by the physician is required prior to use of the device.

Warnings

- Contents are supplied sterile and are intended for single use only. Do not attempt to clean or re-sterilize this product.
- Reuse or reprocessing has not been evaluated and may lead to its failure and subsequent patient illness, infection or other injury.
- Do not use if package is open or damaged and if the expiry date has been exceeded.
- Inspect the package integrity before use.
- Do not continue to use if any of the component are damaged during the procedure.
- The device is intended for continuous use for less than 60 minutes.
- Contact with hard tissue during insertion of the device is considered off-label misuse.
- Use on central nervous system is off-label misuse of the device.
- The device is used for non-vascular procedures; vascular dilation is considered off-label misuse.
- It is deliberate misuse for the patient to remove dilator or any associated accessories.

Precautions

- The product must only be used by a qualified, licensed physician familiar with the percutaneous drainage procedure, with the possible side effects, typical findings, limitations, indications and contraindications of percutaneous drainage procedures.
- Risk is minimal with no long-term effect if device is properly used by a qualified physician.
- During insertion contact with bone, cartilage and scar tissue can damage tip.
- The Screw Dilator can be used to implant a catheter for patients who are ambulatory.

Note: These instructions are not meant to define or suggest any medical or surgical technique. The individual physician is responsible for the proper procedure and techniques to be used with this device.

Directions for Use

- Using the Seldinger Technique, establish access to the tissue and place the guidewire ($\leq 0.038"$ and not included in the package) within the access site.
- Using the Screw Dilator, gradually dilate the tissue to one or two French sizes above the size of the drainage catheter. This will ease the introduction of the drainage catheter.
- Insert the catheter over the wire.
- Remove the guidewire and Screw Dilator leaving the catheter in place.

Disposal

After use, this product may be a potential biohazard. Handle and dispose in accordance with hospitals policies and procedures for waste and biohazard materials.

Storage

Store in controlled room temperature.

NOTE: In the event a serious incident related to this device occurs, the event should be reported to Argon Medical at quality.regulatory@argonmedical.com as well as to the competent health authority where the user/patient resides.

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