


Introducing the redesigned

V·Stick™


Vascular Access Sets

Now with the *UltraSharp* Needle
and *SmoothTrans* Co-axial Introducer



V•Stick™ Vascular Access Sets

UltraSharp Needle - Easy Penetration



**New V•Stick
UltraSharp
Needle**



*Testing performed on polyethylene film using sterile product available on the market

SmoothTrans Co-axial Introducer - Virtually Transitionless



V•Stick™ Vascular Access Sets w/Standard Co-axial Introducer (Box of 10)

Kit Contents: 21 Ga needle (7 cm unless otherwise specified), .018" guidewire, Co-axial introducer (10 cm)

CATALOGUE NO.	DESCRIPTION
391504000E	4F V•Stick™, Standard Introducer, stainless steel guidewire, echogenic
391504300E	4F V•Stick™, Standard Introducer, nitinol guidewire, echogenic
391504400E	4F V•Stick™, Standard Introducer, nitinol guidewire, non-echogenic
391504500E	4F V•Stick™, Standard Introducer, stainless steel guidewire, non-echogenic
391505000E	5F V•Stick™, Standard Introducer, stainless steel guidewire, echogenic
391505300E	5F V•Stick™, Standard Introducer, nitinol guidewire, echogenic
391505400E	5F V•Stick™, Standard Introducer, nitinol guidewire, non-echogenic
391505500E	5F V•Stick™, Standard Introducer, stainless steel guidewire, non-echogenic
390204400E	4F V•Stick™, Standard Introducer, nitinol guidewire, non-echogenic 3.8cm needle

V•Stick™ Vascular Access Sets w/Stiff Co-axial Introducer (Box of 10)

Kit contents: 21 Ga needle (7cm), .018" guidewire, co-axial introducer (10cm)

CATALOGUE NO.	DESCRIPTION
391604000E	4F V•Stick™, Stiff Introducer, stainless steel guidewire, echogenic
391604300E	4F V•Stick™, Stiff Introducer, nitinol guidewire, echogenic
391604400E	4F V•Stick™, Stiff Introducer, nitinol guidewire, non-echogenic
391604500E	4F V•Stick™, Stiff Introducer, stainless steel guidewire, non-echogenic
391605000E	5F V•Stick™, Stiff Introducer, stainless steel guidewire, echogenic
391605300E	5F V•Stick™, Stiff Introducer, nitinol guidewire, echogenic
391605400E	5F V•Stick™, Stiff Introducer, nitinol guidewire, non-echogenic
391605500E	5F V•Stick™, Stiff Introducer, stainless steel guidewire, non-echogenic

Argon Medical Devices

+33 (0)2 99 68 94 95 Customer Service

+33 (0)2 99 23 87 21 Fax

info-europe@argonmedical.com

www.argonmedical.com

2014 Argon Medical Devices, Inc. All Rights Reserved. V•Stick™ is a trademark of Argon Medical Devices, Inc.
96-6300-33/0714A

