

UltraStream

Chronic Dialysis Catheter

THREE LUMENS FOR RELIABLE DIALYSIS CARE




Durable Luer Connectors withstand multiple connections throughout routine therapy


Kink-Resistant Carbothane® Material provides strength, durability and patient comfort

Radiopaque Material enhances visibility under fluoroscopy

Advanced Tip Design maintains zero recirculation during normal flow¹



Dual Arterial Lumens promote continuous flow, even if an arterial lumen is blocked



Smooth Tip Transition facilitates easy insertion over a wire

KIT CONFIGURATIONS

Peel Away Kit	Standard Kit	Exchange Kit
UltraStream Catheter (15.5F)	UltraStream Catheter (15.5F)	UltraStream Catheter (15.5F)
Flexible Stiffener	Flexible Stiffener	Flexible Stiffener
Needleless Injection Site Caps (2)	Needleless Injection Site Caps (2)	Needleless Injection Site Caps (2)
18Ga x 2.75" Introducer Needle	18Ga x 2.75" Introducer Needle	
Vessel Dilator (16F)	Vessel Dilator (16F)	
#11 Safety Scalpel	#11 Safety Scalpel	
Tunneling Device	Tunneling Device	
Adhesive Bandages (2)	Adhesive Bandages (2)	
Safety Foam	Safety Foam	
0.038" x 80cm Guidewire		
Peel Away Sheath (16F)		

ULTRASTREAM ORDERING INFORMATION

CATALOGUE #	DESCRIPTION
32401524	15.5F X 24cm Catheter -- Peel Away Sheath Kit
32401528	15.5F X 28cm Catheter -- Peel Away Sheath Kit
32401532	15.5F X 32cm Catheter -- Peel Away Sheath Kit
32401536	15.5F X 36cm Catheter -- Peel Away Sheath Kit
32401540	15.5F X 40cm Catheter -- Peel Away Sheath Kit
32101524	15.5F X 24cm Catheter -- Standard Kit
32101528	15.5F X 28cm Catheter -- Standard Kit
32101532	15.5F X 32cm Catheter -- Standard Kit
32101536	15.5F X 36cm Catheter -- Standard Kit
32101540	15.5F X 40cm Catheter -- Standard Kit
32001524	15.5F X 24cm Catheter -- Exchange Kit
32001528	15.5F X 28cm Catheter -- Exchange Kit
32001532	15.5F X 32cm Catheter -- Exchange Kit
32001536	15.5F X 36cm Catheter -- Exchange Kit
32001540	15.5F X 40cm Catheter -- Exchange Kit

¹ Data on file.

Carbothane is a registered trademark of the Lubrizol Corporation.

© 2018 Argon Medical Devices, Inc., All Rights Reserved.

No reproduction of any portion of this work shall be made without prior permission of the author.

96-6400-36/0618A . www.argonmedical.eu.com



96-6400-36